

## Product Catalogue





Selecting the right mobility aids to meet your needs is an important decision. Many people use assistive devices to increase or maintain independence. Wheelchairs, rollators, walkers, crutches and walking sticks empower you to go further and fully experience the world around you — on your terms.

For over half a century, Days® has been a recognised global leader in empowering personal mobility. Fully understanding and embracing an individual's mobility needs, Days enhances and advances lives at every turn. When it comes to safety, ease of use, affordable quality, customer support and breadth of product choices, the world turns to Days.

When you select a Days mobility aid, you can expect to experience the difference. Days represents:

## Boundlessness

We share our passion for removing life's boundaries, offering clinicians and individuals a broad range of mobility solutions. With a global reach, Days is always nearby.

## Moving

Days is trust that moves you. We strive to inspire excitement about the possibilities of everything life has to offer. We believe your physical needs should not be a limitation, and we are here to help you propel your life — wherever you wish to go.

## Knowledgeable

Days products and services support you by being built on a solid foundation of quality materials and heritage. With a firm understanding of the mobility industry, we are made to exceed expectations. You can rest assured with Days.

## Compassionate

With every design, Days considers the unique needs and demands of the user. Days knows the products it creates will move both the body and spirit of an individual.



## Go further with Days

Founded in 1968 in Bridgend, Wales, Days has led the industry in mobility for over 50 years. With a broad range of wheelchairs, rollators, walkers, crutches and walking sticks, Days is the brand to which the world turns to advance their lives. Our products are easy to use and ready when you need them while offering affordable quality. With a steadfast focus on safety, Days offers products that customers trust for their mobility needs.

The Days foundation rests on four central pillars: dedication, empowerment, understanding and enabling freedom.

### Dedication

With a steadfast focus on product safety and ease of use, customers and clinicians can trust — and rest assured — knowing our products will take care of their mobility needs. With over a half-century of expertise, Days has a time-honoured reputation that continues to propel customers and their lives forward each day.

### Empowerment

Days goes the distance with a quality, broad range of mobility aids that fit each unique individual for life's ever-changing path. Days helps empower clinicians, distributors and individuals to find the ideal solution for daily freedom of movement and motion.

### Understanding

Days has developed an array of mobility support with every diverse need in mind. With extensive knowledge of the mobility industry, Days applies its empathy and understanding to produce quality, affordable products that support motion.

### Enabling freedom

Navigating through life's twists and turns, Days helps keep the momentum going with a determination to uplift, support and answer each user's needs.



# Mobility Accessories

## Go further with Days

Days offers a broad range of mobility accessories and wet weather gear to ensure that customers have everything they need to keep them safe and organised and dry when they're out and about.

All of our products are made from quality materials with functionality and style in mind.



# Economy Wheelchair Bag

The Days Economy Wheelchair Bag consists of hard-wearing waterproof polyester material with durable loops to fit securely over the handles of a wheelchair and keep it secure whilst moving.

It has a large top zipper main compartment for easy access as well as small feet to ensure stability when standing up.

The bag is specially designed to not interfere with the user pushing the wheelchair and features an additional large pocket on the front for extra storage.

The bag features an extra large storage capacity of 27 L. Size 400 x 450 x 140mm (16 x 18 x 5½").

- Hard-wearing waterproof bag
- Durable loops to secure to wheelchair
- Small feet on the bottom to ensure stability



BLACK



GREY



NAVY

**Code**  
**09118868**  
**7028356**  
**7028357**

**Description**  
Days Economy Wheelchair Bag  
Days Economy Wheelchair Bag  
Days Economy Wheelchair Bag

**Colour**  
● Black  
● Grey  
● Navy

# Deluxe Wheelchair Bag

The stylish Days Deluxe Wheelchair bag features several external pockets large enough to hold clothing or a generous amount of shopping. The inside of the bag is lined and padded to ensure your belongings are safe and secure on the go. The bag consists of hard-wearing, waterproof polyester and features an internal pocket for valuables storage.

The Days Deluxe Wheelchair bag also features two side pockets storing umbrellas or drink bottles. Two durable loops ensure the bag attaches and fits securely to the wheelchair's handles.

The bag also has a rubberised top handle for extra comfort and reflective strips to ensure visibility at night. Rubberised zips allow better grip for users opening and closing the bag.

Maximum safe load 8kg.

Size 400 x 600 x 150mm (15¼ x 23½ x 6").

- Lined and padded bag for extra security
- Padded storage pockets for umbrellas or folding walking sticks
- Waterproof internal pocket



BLACK



GREY



NAVY

**Code**  
**091187707**  
**7028366**  
**7028367**

**Description**  
Deluxe Wheelchair Bag  
Deluxe Wheelchair Bag  
Deluxe Wheelchair Bag

**Colour**  
● Black  
● Grey  
● Navy

# Wheelchair Crutch Bag

A large carry bag specially designed for wheelchairs. The bag features two durable loops that hook over the wheelchair's handles to keep the bag securely in place. The bag has two deep side pouches, accommodating both walking sticks and crutches. The sticks or crutches should always be reduced to their minimum length for safe carriage. The padded main compartment and has a zipped top for secure, waterproof storage. With a 27 L capacity, the bag is large enough to hold clothing or a generous amount of shopping. A reflective strip on the back helps to ensure you remain visible at night. The bag itself is made from hard-wearing, waterproof polyester.

The Days Wheelchair Crutch Bag also features a large external pocket for extra storage and a rubberised top handle for extra strength and comfort when carrying. Max load 7kg. Size 400 x 450 x 140mm (16" x 18" x 6").



- Features two long side pockets to accommodate crutches
- Padded compartment for secure storage
- 27ltr capacity

**Code**  
**091187699**  
**7028360**  
**7028361**

**Description**  
Wheelchair Crutch Bag  
Wheelchair Crutch Bag  
Wheelchair Crutch Bag

**Colour**  
● Black  
● Grey  
● Navy



**BLACK**



**GREY**



**NAVY**

# Economy Scooter Bag

The large waterproof storage bag includes multi-use storage compartments. The bag consists of hard-wearing waterproof polyester material and features a large, zipped main storage compartment and a handy, quick-access pouch on the rear of the bag with a hook-and-loop fastener.

The bag features an extra large storage capacity of 27 L. Size 400 x 450 x 140mm (16 x 18 x 5½").

- Adjustable straps to securely fit over the back of the scooter
- Large zipped main compartment
- 27Ltr capacity



**Code**  
**091187731**  
**7028358**  
**7028359**

**Description**  
Wheelchair Crutch Bag  
Wheelchair Crutch Bag  
Wheelchair Crutch Bag

**Colour**  
● Black  
● Grey  
● Navy



**BLACK**



**GREY**



**NAVY**

# Scooter Crutch Bag

The Days Scooter Crutch Bag is made from hard-wearing waterproof material and features internal padding for added protection. Features a large, zipped main storage compartment with a handy, quick-access pouch on the rear of the bag with a hook-and-loop fastener. Equipped with rubberised handle for extra durability and a reflective strip to ensure safety on the road or streets.

The side pouch holds two crutches or walking sticks. When possible, adjustable crutches or walking sticks should be reduced to minimum height for safe carriage.

Extra large storage capacity of 27 L. Size 400 x 450 x 140mm (16 x 18 x 5½").

- Padded waterproof bag for added protection
- Features two long side pockets to accommodate crutches
- 27Ltr capacity



**Code**  
**091187749**  
**7028362**  
**7028363**

**Description**  
Scooter Crutch Bag  
Scooter Crutch Bag  
Scooter Crutch Bag

**Colour**  
● Black  
● Grey  
● Navy



**BLACK**



**GREY**



**NAVY**

# Deluxe Scooter Bag

This scooter bag is highly padded in all pockets and throughout the main compartment to help ensure the safety of your belongings. The bag itself is made from hard-wearing waterproof material and features three large exterior pockets, two for belongings and a third pocket with a drawstring for a folding walking stick or an umbrella.

The Deluxe Scooter Bag also features an internal clear zipped pouch for wallet storage, reflective strips to increase safety on the streets or roads and a rubberised top handle for extra comfort when transporting the bag. The bag straps can be easily adjusted to sit securely and safely attached to the back of a scooter chair.

Extra large storage capacity of 27 L  
Size 400 x 450 x 140mm (16 x 18 x 5½").

- Lined and padded bag for extra security
- Padded storage pockets for umbrellas or folding walking sticks
- Waterproof internal pocket



**Code**  
**091187753**  
**7028364**  
**7028365**

**Description**  
Deluxe Scooter Bag  
Deluxe Scooter Bag  
Deluxe Scooter Bag

**Colour**  
● Black  
● Grey  
● Navy



**BLACK**



**GREY**



**NAVY**

# Wheelchair Cosy

This fleece blanket is designed for wheelchair users. It zips closed to keep the lower body warm and comfortable. The fleece lining is extraordinarily soft and comfortable to sit on for extended periods of time. Features a hidden inside pocket to store valuables away. The Cosy also has a zip extender so the user can put the zip up without assistance.

## Waterproof and Spill Proof

The outer layer of the Days Wheelchair Cosy is made from waterproof material. This is great for preventing messes from spills. It can also be used to protect the user from light rain.

## Keeping Warm

The Wheelchair Cosy zips closed and has a back panel. It stays in place when in use and fully covers the lower body for a warm and comfortable wheelchair experience.

- Fleece lined blanket that zips closed
- Features a waterproof outer layer
- Zipper has an extra large pull-ring
- Designed with a back panel for extra upper body warmth
- Ideal for keeping the lower body and torso warm



|   |                      |
|---|----------------------|
| A | Length               |
|   | Large 111 cm (43.7") |
|   | Medium 96 cm (37.8") |
|   | Small 87 cm (34.25") |
| B | Width                |
|   | Large 73 cm (28.75") |
|   | Medium 70 cm (27.5") |
|   | Small 67 cm (26.4")  |



FLEECE LINED WITH INSIDE POCKET

| Code      | Size   | Colour     |
|-----------|--------|------------|
| 7028342   | Small  | ● Navy     |
| 7028344   | Small  | ● Burgundy |
| 7028343   | Small  | ● Green    |
| 091119080 | Medium | ● Navy     |
| 7028341   | Medium | ● Burgundy |
| 7028340   | Medium | ● Green    |
| 091169820 | Large  | ● Navy     |
| 7028368   | Large  | ● Burgundy |
| 7028089   | Large  | ● Green    |



BURGUNDY



GREEN



NAVY

# Scooter Cape

This practical, lightweight cape offers all-round protection for both you and your scooter. The cape is made from waterproof nylon with a chest zip and has a transparent panel so you can see the controls and ensure road users can still see your lights.

The cape offers two small openings at the front for the users hands in order to hand change or press for a lift. The openings feature an overlap to prevent water getting into the cape.

The cape comes in medium and small to ensure a perfect fit over the user and the scooter. The cape also features a hood for added protection.

## Waterproof and Road Safe

- Made from waterproof nylon
- Transparent panel for road safety
- Available in small and medium to suit your needs



|          |                                                         |
|----------|---------------------------------------------------------|
| <b>A</b> | Height<br>Medium: 90 cm (35.4")<br>Small: 68 cm (26.8") |
| <b>B</b> | Length<br>Medium: 75 cm (29.5")<br>Small: 68 cm (26.7") |
| <b>C</b> | Width<br>Medium: 45 cm (17.7")<br>Small: 45 cm (17.7")  |



WATERPROOF NYLON  
WITH ZIPPER CHEST



TWO OPENINGS FOR  
EASY HAND ACCESS



BUILT-IN HOOD FOR  
ADDED PROTECTION

Code  
**7028346**  
**091118348**

Size  
Small  
Medium

Colour  
● Navy  
● Navy



NAVY

# Wheelchair Poncho

This poncho is a waterproof sunshade tarp for wheelchair users. Water easily slides off the poncho in the rain and covers most of the body for full protection. The Days Wheelchair Poncho is great for rainy days and is made from waterproof, lightweight nylon.

## For Wheelchair Users

The Days Wheelchair Poncho is shaped to cover the user in a seated position and provides a better fit for wheelchair users. The poncho features two rear vents that go over the wheelchair handles to make it easier to push the wheelchair user.

## Waterproof with Hood

- Waterproof poncho with a zip closure and hood
- Rear vents leave wheelchair handles exposed for caretakers
- Extra length poncho for near full body coverage



|          |                                  |
|----------|----------------------------------|
| <b>A</b> | Front<br>100 cm (39.4")          |
| <b>B</b> | Back<br>68 cm (26.8")            |
| <b>C</b> | Width<br>133 cm (52.4")          |
| <b>D</b> | Hood<br>39 x 25 cm (15.3 x 9.8") |



WATERPROOF NYLON  
WITH ZIPPER CHEST



REAR VENTS  
FOR WHEELCHAIR  
HANDLES



BURGUNDY



GREEN



NAVY

**Code**  
**7028348**  
**7028347**  
**091119049**

**Size**  
Universal  
Universal  
Universal

**Colour**  
● Burgundy  
● Green  
● Navy

# Wheelchair Poncho with Sleeves

Designed with a full length zip, this mobility poncho is quick and easy to put on and can be used for wheelchairs and scooters. The cuffs are adjustable giving a smart finish and preventing the user from getting wet. A snug fitting hood and neckline zip combine to offer maximum protection and comfort. The poncho offers maximum coverage without compromising on manoeuvrability.

## For Wheelchair Users

The Days Wheelchair Poncho is shaped to cover the user in a seated position and provides a better fit for wheelchair users.

## Waterproof and Comfortable


- Made from waterproof material to protect the user and scooter
- Adjustable cuffs on the sleeves for added comfort



|          |                                   |
|----------|-----------------------------------|
| <b>A</b> | Front<br>100 cm (39.4")           |
| <b>B</b> | Back<br>68 cm (26.8")             |
| <b>C</b> | Width<br>133 cm (52.4)            |
| <b>D</b> | Hood<br>39 x 25 cm (15.3 x 9.8")  |
| <b>E</b> | Sleeve<br>72 x 28 cm (28.3 x 11") |



WATERPROOF NYLON  
WITH ZIPPER CHEST



ADJUSTABLE CUFFS  
FOR ADDED COMFORT

Code  
**7028351**  
**7028350**  
**7028349**

Size  
Universal  
Universal  
Universal

Colour  
● Burgundy  
● Green  
● Navy



BURGUNDY



GREEN



NAVY

# Wheelchair Mac

The Days Wheelchair Mac is made from lightweight, waterproof nylon and features an elasticated skirt designed to fit around most wheelchairs. The elasticated band ensures a snug and waterproof fit. The Mac is designed to tuck under the user's legs to keep them dry.

## Waterproof and Lightweight

- Made from lightweight, waterproof nylon
- Features rear vents for wheelchair handles
- The elasticated skirt gives a snug, waterproof fit



|          |                              |
|----------|------------------------------|
| <b>A</b> | Height<br>130 cm (51.2")     |
| <b>B</b> | Hood Length<br>36 cm (14.2") |
| <b>C</b> | Zip Length<br>31 cm (12.2")  |



WATERPROOF NYLON  
WITH ZIPPER CHEST



REAR VENTS FOR  
WHEELCHAIR  
HANDLES



BURGUNDY



GREEN



NAVY

Code  
**7028353**  
**7028352**  
**091119056**

Size  
Universal  
Universal  
Universal

Colour  
● Burgundy  
● Green  
● Navy

# Wheelchair Mac with Sleeves

Seam-sealed for total waterproofing, this Mac offers complete protection for wearer and the wheelchair. Features knitted cuffs on the sleeves to prevent the users arm from getting wet. Vents at the rear of the Mac allow the use of the wheelchairs, pushing handles, and a small pouch at the front provides handy storage. Elasticated all the way around to ensure a snug fit.

## Waterproof and Spill Proof

- Seam-sealed for total waterproofing
- Knitted cuffs for added comfort
- Vents at the rear for access to the wheelchair handle



|          |                                |
|----------|--------------------------------|
| <b>A</b> | Height<br>120 cm (47.25")      |
| <b>B</b> | Hood Length<br>36 cm (14.2")   |
| <b>C</b> | Sleeve Length<br>68 cm (26.8") |
| <b>D</b> | Width<br>48 cm (18.9")         |



**WATERPROOF NYLON WITH ZIPPER CHEST**



**REAR VENTS FOR WHEELCHAIR HANDLES**



**KNITTED CUFFS FOR ADDED COMFORT**

**Code**  
**7028355**  
**7028354**  
**091119064**

**Size**  
Universal  
Universal  
Universal

**Colour**  
● Burgundy  
● Green  
● Navy



**BURGUNDY**



**GREEN**



**NAVY**

# Scooter Storage Cover

The Days Scooter Storage cover is windproof and has an elasticated bottom to ensure it can securely tuck under the scooter when parked.


The cover is also waterproof to prevent rainwater from seeping through the material when storing the scooter. The storage cover can easily be wiped clean and is suitable for scooters, wheelchairs and power chairs.

Its universal design allows it to fit most scooters. Please check the dimensions of your scooter and the cover before purchasing. The Days Scooter Cover is easy to store and folds up for easy portability. To wash, rinse with water and wipe dry after use.

- **Elasticated bottom to ensure a good fit**
- **Waterproof design for outdoor storage**
- **Available in small, medium and large to suit your needs**



|          |                                                                                 |
|----------|---------------------------------------------------------------------------------|
| <b>A</b> | Height<br>Large 110 cm (43.3")<br>Medium 104 cm (40.9")<br>Small 87 cm (32.25") |
| <b>B</b> | Length<br>Large 154 cm (60.6")<br>Medium 147 cm (57.9")<br>Small 115 cm (45.3") |
| <b>C</b> | Width<br>Large 71 cm (27.9")<br>Medium 70 cm (27.5")<br>Small 59 cm (23.2")     |



**ELASTICATED BOT-  
TOM ENSURES  
A GOOD FIT**

| Code             | Size   | Colour |
|------------------|--------|--------|
| <b>7028345</b>   | Small  | ● Navy |
| <b>091118314</b> | Medium | ● Navy |
| <b>7028730</b>   | Large  | ● Navy |



**NAVY**



# Wheelchairs

## You're in the driver's seat

The World Health Organization recognizes that wheelchairs are essential for mobility, dignity and general wellbeing<sup>1</sup>. Whether self- or attendant-propelled, your high-quality chair should enhance your mobility and safety while being easy to use and affordable. With Days, you will also have access to customer support when you need it most.

### All Days Wheelchairs include the following features:

- Puncture-proof tyres to help ensure trouble-free mobility, low maintenance and a smooth ride over most surfaces
- Moulded armrests for extra comfort
- Folding for easy storage and transportation
- Wipe-clean padded seat fabric



References: 1. Strengthening access to quality wheelchair services. World Health Organization. <https://www.who.int/teams/health-product-policy-and-standards/assistive-and-medical-technology/assistive-technology/wheelchair-services>. Accessed May 27, 2022.

# Escape Lite Attendant Wheelchair

## Ease of transportation and storage

Escape Lite Wheelchairs are ultra-lightweight wheelchairs suitable for indoors and outdoors. Experience peace of mind and comfortability with puncture-proof tyres, a full-length armrest, padded upholstery and detachable, swing-away footrests. You will have the flexibility to get around no matter where you go.



## Key features

- Lightweight, foldable frame with foldable backrest
- Swing-away and detachable footplates
- Padded, flame-retardant upholstery
- Includes a cable braking system
- Seatbelt included

## Specifications

|                                |                                                                     |
|--------------------------------|---------------------------------------------------------------------|
| Seat width                     | Narrow: 41 cm (16"),<br>Standard: 46 cm (18"),<br>Wide: 51 cm (18") |
| Seat depth                     | 41 cm (16")                                                         |
| Seat height (front)            | 51.5 cm (20¼")                                                      |
| Seat height (rear)             | 49.5 cm (19½")                                                      |
| Width between arms             | Narrow: 38 cm (15")<br>Standard: 41 cm (16")<br>Wide: 48 cm (19")   |
| Backrest height above seat     | 40 cm (15¾")                                                        |
| Footrest length                | 39-49 cm (15-19")                                                   |
| Armrest height above seat      | 23 cm (9")                                                          |
| External width (open)          | Narrow: 52 cm (20½")<br>Standard: 53 cm (21")<br>Wide: 56 cm (22")  |
| External width (closed)        | 23 cm (9")                                                          |
| External depth (w/o footrests) | 75 cm (29½")                                                        |
| External depth (w/ footrests)  | 99.5 cm (39")                                                       |
| Chair weight                   | 12 kg (26 lb)                                                       |
| Chair weight (w/o footrests)   | 10.5 kg (23 lb)                                                     |



| Code             | Description           | Width    | Colour         |
|------------------|-----------------------|----------|----------------|
| <b>091547561</b> | Escape Lite Attendant | Narrow   | ● Purple       |
| <b>091555481</b> | Escape Lite Attendant | Narrow   | ● Racing Green |
| <b>091556901</b> | Escape Lite Attendant | Narrow   | ● Ruby Red     |
| <b>091311521</b> | Escape Lite Attendant | Narrow   | ● Silver Blue  |
| <b>091547579</b> | Escape Lite Attendant | Standard | ● Purple       |
| <b>091555473</b> | Escape Lite Attendant | Standard | ● Racing Green |
| <b>091556893</b> | Escape Lite Attendant | Standard | ● Ruby Red     |
| <b>091171727</b> | Escape Lite Attendant | Standard | ● Silver Blue  |
| <b>091556604</b> | Escape Lite Attendant | Wide     | ● Purple       |
| <b>091556612</b> | Escape Lite Attendant | Wide     | ● Racing Green |
| <b>091556596</b> | Escape Lite Attendant | Wide     | ● Ruby Red     |
| <b>091556588</b> | Escape Lite Attendant | Wide     | ● Silver Blue  |



NARROW,  
STANDARD



WIDE



# Escape Lite Self-Propelled Wheelchair

## Ease of transportation and storage



Escape Lite Wheelchairs are ultra-lightweight wheelchairs suitable for indoors and outdoors. Experience peace of mind and comfortability with puncture-proof tyres, a full-length armrest, padded upholstery and detachable, swing-away footrests. You will have the flexibility to get around no matter where you go.

## Key features



- Lightweight with foldable aluminium frame
- Swing-away and detachable footplates
- Padded, flame-retardant upholstery
- Seatbelt included
- Quick-release MAG wheels
- Attendant and user brakes

## Specifications

|                                |                                                                     |
|--------------------------------|---------------------------------------------------------------------|
| Seat width                     | Narrow: 41 cm (16")<br>Standard: 46 cm (18")<br>Wide: 51 cm (20")   |
| Seat depth                     | 41 cm (16")                                                         |
| Seat height (front)            | 51.5 cm (20¾")                                                      |
| Seat height (rear)             | 49.5 cm (19½")                                                      |
| Backrest height above seat     | 40 cm (15¾")                                                        |
| Armrest height above seat      | 23 cm (9")                                                          |
| Footrest length                | 39-49 cm (15-19")                                                   |
| External width (open)          | Narrow: 62 cm (24")<br>Standard: 64.5 cm (25")<br>Wide: 69 cm (27") |
| External width (closed)        | 31 cm (12¼")                                                        |
| External depth (w/o footrests) | 80 cm (31½")                                                        |
| External length (w/ footrests) | 106 cm (41½")                                                       |
| External height                | 94 cm                                                               |
| Chair weight                   | 15.5 kg (34 lb)                                                     |
| Chair weight (w/o footrests)   | 14 kg (31 lb)                                                       |



| Code             | Description                | Width    | Colour         |
|------------------|----------------------------|----------|----------------|
| <b>091566314</b> | Escape Lite Self-Propelled | Narrow   | ● Purple       |
| <b>091566298</b> | Escape Lite Self-Propelled | Narrow   | ● Ruby Red     |
| <b>091566280</b> | Escape Lite Self-Propelled | Narrow   | ● Silver Blue  |
| <b>091566272</b> | Escape Lite Self-Propelled | Standard | ● Purple       |
| <b>091566264</b> | Escape Lite Self-Propelled | Standard | ● Racing Green |
| <b>091566256</b> | Escape Lite Self-Propelled | Standard | ● Ruby Red     |
| <b>091566249</b> | Escape Lite Self-Propelled | Standard | ● Silver Blue  |
| <b>7500092</b>   | Escape Lite Self-Propelled | Wide     | ● Purple       |
| <b>7500091</b>   | Escape Lite Self-Propelled | Wide     | ● Silver Blue  |



# Swift Aluminium Wheelchair

**SELF OR ATTENDANT  
PROPELLED**

## Versatile and sturdy

Swift Wheelchairs are versatile, easy-to-use chairs. Lightweight with a robust and high-quality aluminium frame, Swift can be quickly disassembled without tools, making transportation and storage simple. Swift is compact, manoeuvrability and ideal for navigating shops, restaurants and tight spaces.

## Key features

- Designed for indoor or outdoor occasional use
- Half-folding back for compact storage and convenient reaching
- Attendant-propelled version includes a cable braking system






## Specifications

| Seat width                     | 41 cm (16"), 46 cm (18") |                            |
|--------------------------------|--------------------------|----------------------------|
| Seat depth                     | 42 cm (16.5")            |                            |
| Seat height (front)            | 51 cm (20")              |                            |
| Seat height (rear)             | 48.5 cm (19")            |                            |
| Backrest height above seat     | 40 cm (15¾")             |                            |
| Footrest length                | 44-54 cm (17¼-21¼")      |                            |
| Armrest height above seat      | 23 cm (9")               |                            |
|                                | Self-propelled           | Attendant-propelled        |
| External width (open)          | 62 cm (24"), 67 cm (26") | 57 cm (22½"), 62 cm (24½") |
| External width (closed)        | 28 cm (11")              | 24 cm (9½")                |
| External depth (w/o footrests) | 83 cm (32¾")             | 75 cm (29½")               |
| External depth (w/ footrests)  | 113 cm (44½")            | 107 cm (42")               |
| Chair weight                   | 12.5 kg (27.5 lb)        | 11.5 kg (25.3 lb)          |
| Chair weight (w/o footrests)   | 11.2 kg (24.7 lb)        | 10.2 kg (22.4 lb)          |



| Code             | Description               | Width      | Colour        |
|------------------|---------------------------|------------|---------------|
| <b>091186618</b> | Swift Attendant Propelled | 41cm (16") | ● Cool Silver |
| <b>091186626</b> | Swift Attendant Propelled | 46cm (18") | ● Cool Silver |
| <b>091555515</b> | Swift Attendant Propelled | 46cm (18") | ● Turquoise   |
| <b>091555523</b> | Swift Attendant Propelled | 46cm (18") | ● Burgundy    |
| <b>091186683</b> | Swift Self Propelled      | 46cm (18") | ● Turquoise   |
| <b>091555556</b> | Swift Self Propelled      | 46cm (18") | ● Burgundy    |



# Whirl Steel Wheelchair

**SELF OR ATTENDANT  
PROPELLED**

## Ideal for part- or full-time use

Whirl Wheelchairs are simple yet versatile and ideal for part or full-time use. Whirl has a detachable, half-folding backrest for compact storage, and its narrow frame and compact dimensions ensure excellent manoeuvrability indoors and confined spaces, and tight spaces.



## Key features

- Removable desk armrests to enable closer access to tables and worktops
- Self-propelled version features removable quick-release rear wheels
- Passes ISO 7176 part 19






## Specifications

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| Seat width                 | 41 cm (16"), 43 cm (17"), 45 cm (17¾"), 48 cm (19") |
| Seat depth                 | 42 cm (16½")                                        |
| Seat height (front)        | 51 cm (20")                                         |
| Seat height (rear)         | 48 cm (19")                                         |
| Backrest height above seat | 43 cm (17")                                         |
| Footrest length            | 43.5-50 cm (17¼"-19½")                              |
| Armrest height above seat  | 23 cm (9")                                          |

|                                | Self-propelled           | Attendant-propelled     |
|--------------------------------|--------------------------|-------------------------|
| External width (open)          | Seat width +20 cm (7¾"), | Seat width +16 cm (6¾") |
| External width (closed)        | 29 cm (11½")             | 25 cm (9¾")             |
| External depth (w/o footrests) | 80 cm (31½")             | 70 cm (27½")            |
| External depth (w/ footrests)  | 105 cm (41¾")            | 95 cm (37½")            |
| Chair weight                   | 17 kg (37.4 lb)          | 15.4 kg (33.9 lb)       |
| Chair weight (w/o footrests)   | 15 kg (33 lb)            | 13.4 kg (29.5 lb)       |



| Code             | Description               | Width       | Colour   |
|------------------|---------------------------|-------------|----------|
| <b>091440346</b> | Whirl Attendant Propelled | 41cm (16")  | ● Silver |
| <b>091440361</b> | Whirl Attendant Propelled | 43cm (17")  | ● Silver |
| <b>091440379</b> | Whirl Attendant Propelled | 45cm (17¾") | ● Silver |
| <b>091440395</b> | Whirl Attendant Propelled | 48cm (19")  | ● Silver |
| <b>091440270</b> | Whirl Self Propelled      | 41cm (16")  | ● Silver |
| <b>091440288</b> | Whirl Self Propelled      | 43cm (17")  | ● Silver |
| <b>091440296</b> | Whirl Self Propelled      | 45cm (17½") | ● Silver |
| <b>091440304</b> | Whirl Self Propelled      | 45cm (19")  | ● Silver |
| <b>091467729</b> | Optional Leg Rest - Left  |             |          |
| <b>091467760</b> | Optional Leg Rest - Right |             |          |



# Chrome Plated Steel Wheelchair

**SELF OR ATTENDANT  
PROPELLED**

## Robust and resilient

Chrome-plated Steel Wheelchairs, designed for resiliency, are manufactured from heavy gauge steel and have an attractive chrome finish. The embossed vinyl upholstery is hard-wearing and easy to wipe clean.

## Key features

- Stainless steel skirt guards
- Chrome finish
- Swing-away removable arms
- Optional elevating leg rests







## Specifications

| Seat width                     | Narrow: 41 cm (16"),<br>Standard: 46 cm (18") |                            |
|--------------------------------|-----------------------------------------------|----------------------------|
| Seat depth                     | 41 cm (16")                                   |                            |
| Seat height (front)            | 51 cm (20")                                   |                            |
| Seat height (rear)             | 47 cm (18½")                                  |                            |
| Backrest height above seat     | 41 cm (16")                                   |                            |
| Footrest length                | 38-53 cm (15-21")                             |                            |
| Armrest height above seat      | 25 cm (9¾")                                   |                            |
|                                | Self-propelled                                | Attendant-propelled        |
| External width (open)          | 62 cm (21½"), 67 cm (26¾")                    | 58 cm (22¾"), 63 cm (24¾") |
| External width (closed)        | 30 cm (11¾")                                  | 26 cm (10¾")               |
| External depth (w/o footrests) | 80 cm (31½")                                  | 73 cm (28¾")               |
| External depth (w footrests)   | 105 cm (41¾")                                 | 95 cm (37½")               |
| Chair weight                   | 20 kg (44 lb)                                 | 17 kg (37.4 lb)            |
| Chair weight (w/o footrests)   | 18 kg (39.6 lb)                               | 15 kg (33 lb)              |



| Code             | Description                                            | Width                | Back         | Colour   |
|------------------|--------------------------------------------------------|----------------------|--------------|----------|
| <b>091439330</b> | Self-Propelled Wheelchair                              | Standard 46 cm (18") | Fixed Back   | ● Chrome |
| <b>091439132</b> | Self-Propelled Wheelchair                              | Standard 46 cm (18") | Folding Back | ● Chrome |
| <b>091439223</b> | Self-Propelled Wheelchair                              | Narrow 41 cm (16")   | Fixed Back   | ● Chrome |
| <b>091439173</b> | Self-Propelled Wheelchair                              | Narrow 41 cm (16")   | Folding Back | ● Chrome |
| <b>091439116</b> | Attendant-Propelled Wheelchair                         | Standard 46 cm (18") | Fixed Back   | ● Chrome |
| <b>091356708</b> | Attendant-Propelled Wheelchair                         | Standard 46 cm (18") | Folding Back | ● Chrome |
| <b>091439397</b> | Attendant-Propelled Wheelchair                         | Narrow 41 cm (16")   | Fixed Back   | ● Chrome |
| <b>091439959</b> | Days Self-Propelled Wheelchair                         | Narrow 41 cm (16")   | Folding Back | ● Chrome |
| <b>091444439</b> | Chrome-plated elevating leg rests to suit heavy chairs |                      |              |          |



# Heavy-duty Steel Wheelchair

SELF OR ATTENDANT  
PROPELLED

## Designed for strength and durability

Heavy Duty Steel Wheelchairs are manufactured from heavy gauge steel with double crossbar frames and reinforced double-thickness upholstery. Height-adjustable, swing-away and removable footrests come with a calf strap and rubberised footplates for safety.

## Key features






- Stainless steel skirt guards
- Chrome finish
- Optional elevating leg rests

## Specifications

| Seat width                     | Wide: 51 cm (20")<br>Extra wide: 56 cm (22")   |                     |
|--------------------------------|------------------------------------------------|---------------------|
| Seat depth                     | 41 cm (16")                                    |                     |
| Seat height (front)            | 52 cm (20½")                                   |                     |
| Seat height (rear)             | 47 cm (18½")                                   |                     |
| Backrest height above seat     | 41 cm (16")                                    |                     |
| Footrest length                | 40-53 cm (16-21")                              |                     |
| Armrest height above seat      | 25 cm (9¾")                                    |                     |
| External width (open)          | Wide: 68 cm (26¾")<br>Extra wide: 73 cm (28¾") |                     |
| External width (closed)        | 30 cm (11¾")                                   |                     |
|                                | Self-propelled                                 | Attendant-propelled |
| External depth (w/o footrests) | 82 cm (32¼")                                   | 75 cm (29½")        |
| External depth (w/ footrests)  | 107 cm (42")                                   | 100 cm (39½")       |
| Chair weight                   | 25 kg (55 lb)                                  | 23 kg (50.6 lb)     |
| Chair weight (w/o footrests)   | 23 kg (50.6 lb)                                | 21 kg (46.2 lb)     |



| Code             | Description                                            | Width                  | Back         | Colour   |
|------------------|--------------------------------------------------------|------------------------|--------------|----------|
| <b>091356674</b> | Self-Propelled Wheelchair                              | Wide 51 cm (20")       | Fixed Back   | ● Chrome |
| <b>091439181</b> | Self-Propelled Wheelchair                              | Wide 51 cm (20")       | Folding Back | ● Chrome |
| <b>091439280</b> | Self-Propelled Wheelchair                              | Extra-wide 56 cm (22") | Fixed Back   | ● Chrome |
| <b>091227412</b> | Attendant-Propelled Wheelchair                         | Wide 51 cm (20")       | Folding Back | ● Chrome |
| <b>091227420</b> | Attendant-Propelled Wheelchair                         | Extra-wide 56 cm (22") | Fixed Back   | ● Chrome |
| <b>091444439</b> | Chrome-plated elevating leg rests to suit heavy chairs |                        |              |          |



# Fully Reclining Steel Wheelchair

## Fully reclining for increased comfort

Fully Reclining Steel Wheelchairs feature a hydraulic recline mechanism to support your weight so that you can recline with infinite adjustment from vertical to horizontal. Removable head support ensures comfort at every angle. The variation in angle provides the user with a change in orientation, enabling redistribution of pressure to increase comfort and sitting tolerance and reduce the risk of pressure sores. With wipe-clean seat fabric, these chairs are durable with minimum ongoing maintenance.

## Key features

- Self-propelling wheels
- Chrome-plated steel frame
- Stainless steel skirt guards
- Detachable armrests and elevating leg rests

## Specifications

|                                |                                              |
|--------------------------------|----------------------------------------------|
| Seat width                     | Narrow: 41 cm (16")<br>Standard: 46 cm (18") |
| Seat depth                     | 43 cm (17")                                  |
| Seat height (front)            | 51 cm (20")                                  |
| Seat height (rear)             | 47 cm (18½")                                 |
| Backrest height above seat     | 54 cm (22¼")                                 |
| Footrest length                | 40-53 cm (15¾-21")                           |
| Armrest height above seat      | 25 cm (9¾")                                  |
| External width (open)          | Narrow: 61 cm (24")<br>Standard: 66 cm (26") |
| External width (closed)        | 27 cm (10½")                                 |
| External depth (w/o footrests) | 91 cm (35¾")                                 |
| External depth (w/ footrests)  | 125 cm (49¾")                                |
| Chair weight                   | 28 kg (61.6 lb)                              |
| Chair weight (w/o footrests)   | 24 kg (52.8 lb)                              |



COMPACT FOLD  
FOR PORTABILITY

**Code**  
**091227644**  
**091439868**

**Description**  
Self-Propelled Wheelchair  
Self-Propelled Wheelchair

**Width**  
Standard 46 cm (18")  
Narrow 41 cm (16")

**Colour**  
● Chrome  
● Chrome

115 kg

253½ lb

Maximum user weight

1 YEAR

Product warranty

2 YEAR

Main frame warranty

Self propelled


# Rollators

## Your guide to choosing the perfect rollator

Choosing a rollator with a frame at the correct height is essential to ensure wellbeing. Other factors to consider include where you will primarily use your rollator, the amount of support you will need, the number of wheels (three or four), portability and storage.

### When taking measurements, consider the following:

- If the frame is too high, you may find it difficult to straighten your elbows and increase weight distribution through your arms.
- Your frame is too low if you are bent over with poor posture.
- Generally, the height of the handgrips should be at the level of the wrist bone when your elbows are slightly bent, at an angle of about 15° flexion.
- When measuring the height of your walking frame, wear appropriate and supportive footwear.
- If your rollator includes a seat, consider the seat height. When seated comfortably, your feet should be flat on the ground.



# Lightweight 100 Series Rollator


## Keeping you on the move

The 100 Series Rollator is equipped with aluminium frames, combining lightweight design and robust engineering quality with a maximum user weight of 165 kg (26 st or 364¾ lb).

## Key features

- Available in four sizes to accommodate a range of heights and builds
- Lockable loop brakes for security and convenience
- Adjustable brake tension to suit your needs
- Folds easily for storage and transport
- Ergonomic handgrips and adjustable height handles with five settings at 2.5-cm intervals

102-103 series



|                                  | 150-175 cm<br>(5'0"-5'9")                   | 155-180 cm<br>(5'1"-5'11")              | 168-188 cm<br>(5'5"-6'2")             | 175-200 cm<br>(5'9"-6'5")                   |
|----------------------------------|---------------------------------------------|-----------------------------------------|---------------------------------------|---------------------------------------------|
|                                  | X-Small                                     | Small                                   | Medium                                | Large                                       |
| A                                | Handle height<br>77.5-87.5 cm<br>(30½-34½") | Handle height<br>78-90 cm<br>(30¾-35½") | Handle height<br>84-94 cm<br>(33-37") | Handle height<br>87.5-99.5 cm<br>(34½-39¼") |
| B                                | Seat height<br>50 cm (19¾")                 | Seat height<br>53.5 cm (21")            | Seat height<br>56.5 cm (22¼")         | Seat height<br>62 cm (24¼")                 |
| C                                | Backrest height<br>75 cm (29½")             | Backrest height<br>80 cm (31½")         | Backrest height<br>80 cm (31½")       | Backrest height<br>87 cm (34¼")             |
| D                                | Wheel diameter<br>15 cm (6")                | Wheel diameter<br>20 cm (8")            | Wheel diameter<br>20 cm (8")          | Wheel diameter<br>20 cm (8")                |
| Seat width 35 cm (13¾")          |                                             |                                         | Seat depth 32 cm (12¾")               |                                             |
| Product weight<br>6 kg (13.2 lb) |                                             | Product weight<br>7 kg (15.4 lb)        | Product weight<br>7.2 kg (15.8 lb)    | Product weight<br>7.3 kg (16 lb)            |


| Code      | Description                      | Colour       |
|-----------|----------------------------------|--------------|
| 091553023 | 102 X-Small Lightweight Rollator | Blue         |
| 091553072 | 102 X-Small Lightweight Rollator | Purple       |
| 091553080 | 102 X-Small Lightweight Rollator | Quartz       |
| 091356302 | 102 X-Small Lightweight Rollator | Racing Green |
| 091553049 | 102 X-Small Lightweight Rollator | Ruby Red     |

|           |                                |              |
|-----------|--------------------------------|--------------|
| 091553114 | 103 Small Lightweight Rollator | Blue         |
| 091553130 | 103 Small Lightweight Rollator | Purple       |
| 091553098 | 103 Small Lightweight Rollator | Quartz       |
| 091553148 | 103 Small Lightweight Rollator | Racing Green |
| 091553106 | 103 Small Lightweight Rollator | Ruby Red     |

|           |                                 |               |
|-----------|---------------------------------|---------------|
| 091356336 | 105 Medium Lightweight Rollator | Blue          |
| 091546563 | 105 Medium Lightweight Rollator | Purple        |
| 091553163 | 105 Medium Lightweight Rollator | Quartz        |
| 091553171 | 105 Medium Lightweight Rollator | Racing Green  |
| 091356328 | 105 Medium Lightweight Rollator | Ruby Red      |
| 091553155 | 105 Medium Lightweight Rollator | Russet Orange |

|           |                                |               |
|-----------|--------------------------------|---------------|
| 081561661 | 106 Large Lightweight Rollator | Blue          |
| 091553213 | 106 Large Lightweight Rollator | Purple        |
| 091553197 | 106 Large Lightweight Rollator | Quartz        |
| 091553221 | 106 Large Lightweight Rollator | Racing Green  |
| 091553189 | 106 Large Lightweight Rollator | Ruby Red      |
| 091553205 | 106 Large Lightweight Rollator | Russet Orange |

105-106 series




# Four-wheeled Rollator

## Quality and value

Four-wheeled Rollators combine superb quality with fantastic value and affordability. Each rollator features height-adjustable handles, safe-stop, lockable brakes, large swivel wheels and a lightweight aluminium frame.

## Key features

- Affordable and sturdy
- Padded seat and backrest for comfort
- Includes wire basket that hangs under the seat or hooks onto the front of the frame
- Folds easily for transport and storage



**A** Handle height  
78-90 cm (30-35½")

**B** Overall depth  
68 cm (26¾")

**C** Seat height  
26.5 cm (10½")

**D** Backrest height  
26.5 cm (10½")

**E** Wheel diameter  
20 cm (8")

Seat width  
33 cm (13")

Seat depth  
31.5 cm (12½")

Product weight  
6.8 kg (15 lb)



**Code**  
**091169945**  
**091169952**

**Description**  
Four-Wheeled Rollator  
Four-Wheeled Rollator

**Colour**  
● Silver  
● Quartz

165 kg  
26 st 363¼ lb  
Maximum user weight


# Deluxe Lightweight Rollator

## Modern and compact

Deluxe Lightweight Rollators are sleek and stylish. This rollator can be folded into a 28-cm standing position to fit through tight spaces, with an additional option to fold it further into a compact size (28 cm high x 28 cm wide x 39 cm deep).

## Key features

- Lockable loop brakes for security and convenience
- Kerb climbers: Step down on the kerb climber while gently pulling back on the handles to raise the front wheels and easily manoeuvre over obstacles
- Walking stick holder
- Built-in rest seat with canvas backrest



|                                 |                                                                 |
|---------------------------------|-----------------------------------------------------------------|
| <b>A</b>                        | Handle height<br>81-93 cm (32¾-36½")                            |
| <b>B</b>                        | Overall depth<br>68 cm (26¾")                                   |
| <b>C</b>                        | Seat height<br>55 cm (21¾")                                     |
| <b>D</b>                        | Backrest height<br>20.5-33 cm (8-13")                           |
| <b>E</b>                        | Wheel diameter<br>Front - 19.5 cm (7⅞")<br>Rear - 16.5 cm (6½") |
| Seat width<br>46 cm (18")       |                                                                 |
| Seat depth<br>20 cm (7¾")       |                                                                 |
| Product weight<br>7 kg (15½ lb) |                                                                 |



FOLDS INTO A COMPACT SIZE, 66 CM HIGH, 28 CM WIDE, 39 CM DEEP



FOLDS INTO A STANDING POSITION, 28 CM WIDE



LOCKABLE LOOP BRAKES

Code  
**091439314**

Description  
Deluxe Lightweight Rollator

Colour  
● Silver

130 kg  
28½ st 286¾ lb  
Maximum user weight


# Essentials 252 Rollator

## Ideal for everyday use

Essentials 252 Rollators are entry-level, affordable rollators with lightweight frames, more compact seats and lower user weights than the 100 Series, making these rollators ideal for smaller users.

## Key features

- Lockable loop brakes for security and convenience
- Padded seat and backrest
- Ergonomic handgrips for optimum comfort
- Folds easily for transport and storage



|          |                                           |
|----------|-------------------------------------------|
| <b>A</b> | Handle height<br>879.5-90.5 cm (31¼-35½") |
| <b>B</b> | Overall depth<br>68.5 cm (27")            |
| <b>C</b> | Seat height<br>53 cm (20¾")               |
| <b>D</b> | Backrest height<br>30.5 cm (12")          |
| <b>E</b> | Wheel diameter<br>20 cm (8")              |
|          | Seat width<br>30 cm (11¾")                |
|          | Seat depth<br>30.5 cm (12")               |
|          | Product weight<br>6.42 kg (14¼ lb)        |



**Code**  
**091433713**

**Description**  
Deluxe Lightweight Rollator

**Colour**  
● Silver

100 kg  
153¼ st 220½ lb  
Maximum user weight


# Steel Four-wheeled Rollator

## Durable for everyday use

The Steel Four-wheeled Rollator is a durable, height-adjustable walker manufactured from steel tubing with four solid tyre wheels, ideal for anyone who needs a mobility aid. The rest seat and backrest are padded, and the rollator includes locking brakes and a convenient storage bag. It can be folded flat for storage or transportation.

## Key features

- Steel tubing frame with solid wheels
- Lockable loop brakes
- Padded seat and backrest
- Built-in rest seat with canvas backrest



|          |                                                                |
|----------|----------------------------------------------------------------|
| <b>A</b> | Handle height<br>81.2-94.4 cm (32-37")                         |
| <b>B</b> | Overall depth<br>70 cm (27½")                                  |
| <b>C</b> | Seat height<br>Front — 53.9 cm (21¼")<br>Rear — 54.4 cm (21½") |
| <b>D</b> | Backrest height<br>30.5 cm (12")                               |
| <b>E</b> | Wheel diameter<br>20 cm (8")                                   |
|          | Seat width<br>34.5 cm (13½")                                   |
|          | Seat depth<br>29.5 cm (11½")                                   |
|          | Product weight<br>8.2 kg (18 lb)                               |

Code  
**091560333**

Description  
Steel Four-Wheeled Rollator

Colour  
● Red

160 kg  
26 st 352¼ lb  
Maximum user weight


# Heavy-duty Rollator

## Built for strength

Heavy-duty Rollators are manufactured with gauge steel tubing and offer wider frames and rest seats for users who need more space.

## Key features

- Lockable loop brakes for security and convenience
- Under-seat basket for functional storage space
- High user weights
- Folds easily for transport and storage



|          |                                     |
|----------|-------------------------------------|
| <b>A</b> | Handle height<br>83-94 cm (32¾-37") |
| <b>B</b> | Overall depth<br>66 cm (26")        |
| <b>C</b> | Overall width<br>67 cm (26½")       |
| <b>D</b> | Seat height<br>55 cm (21½")         |
| <b>E</b> | Wheel diameter<br>20 cm (8")        |
|          | Seat width<br>49 cm (19¼")          |
|          | Seat depth<br>35 cm (13¾")          |
|          | Product weight<br>13 kg (28½ lb)    |



**Code**  
**091439090**

**Description**  
Heavy-Duty Rollator

**Colour**  
● Blue

178 kg  
28 st 392¼ lb  
Maximum  
user weight


# 111 Bariatric Rollator

## Built for strength

The 111 Bariatric Rollator, constructed of high-quality steel, was designed to help users maintain balance and prevent falls while navigating indoor and outdoor surfaces. The rollator is height-adjustable with a padded backrest to prevent the user from tilting backward when seated. Complete with a storage basket for personal belongings, this rollator folds for easy storage.

## Key features

- High-quality steel frame with height-adjustable handgrips
- Lockable loop brakes with adjustable tension
- Extra-wide padded seat and backrest



**A** Handle height  
83-96 1/2 cm (33-38")

**B** Overall depth  
71 cm (28")

**C** Overall width  
61 cm (24")



**D** Seat height  
55 cm (21 1/2")

**E** Wheel diameter  
20 cm (8")

Seat width  
49 cm (19 3/4")

Seat depth  
35 cm (13 3/4")

Product weight  
13 kg (28 1/2 lb)



Code  
**7500348**

Description  
111 Bariatric Rollator

Colour  
● Red

317 1/2 kg  
50 st 700 lb  
Maximum user weight


# 2-in-1 Rollator & Transit Chair

## Comfort when you need it

The 2-in-1 Rollator and Transit Chair offers the comfort you need for various uses. Each rollator includes a padded seat and backrest, an under-seat bag for storage and flip-up footrests to provide extra support.

## Key features

- Suitable for use as a rollator or transit chair
- Lockable loop brakes for security and convenience
- Folds easily for transport and storage



|          |                                         |
|----------|-----------------------------------------|
| <b>A</b> | Handle height<br>80-93 cm (31¾-36½")    |
| <b>B</b> | Pushing handle height<br>90.5 cm (35½") |
| <b>C</b> | Overall depth<br>74 cm (29¼")           |
| <b>D</b> | Overall width<br>58.5 cm (23")          |
| <b>E</b> | Seat height<br>56 cm (22")              |
| <b>F</b> | Wheel diameter<br>10 cm (8")            |
|          | Seat width<br>46 cm (18¼")              |
|          | Seat depth<br>33 cm (13")               |
|          | Product weight<br>8.8 kg (19¼ lb)       |



**Code**  
**091207299**  
**091207307**

**Description**  
2-in-1 Rollator and Transit Chair  
2-in-1 Rollator and Transit Chair

**Colour**  
● Silver  
● Ruby Red

160 kg  
26 st 352¾ lb  
Maximum user weight

1 YEAR  
Product warranty

2 YEAR  
Main frame warranty


# Breeze Indoor Rollator

## A smooth journey from room to room

Breeze Indoor Rollators offer a stylish alternative to a household trolley. Manufactured from painted aluminium tubing with four solid tyre wheels, each rollator includes a storage bag, locking brakes and height-adjustable handles.

### Key features

- One-handed braking with brake lock
- Easily adjustable handle height
- Durable plastic tray and carry bag with zip pocket
- Fold button for easy storage



**A** Handle height  
83.5-94.5 cm (32¾-37¼")

**B** Overall depth  
66 cm (30")

**C** Overall width  
54 cm (21¼")

**D** Wheel diameter  
17.75 cm (7")

**E** Tray depth  
17.75 cm (7")

Tray width  
49 cm (19¼")

Tray height  
68 cm (26¾")


Tray max load  
5 kg (11 lb)

Bag max load  
5 kg (11 lb)

Product weight  
5.5 kg (12 lb)



**BRAKE LOCK**



**EASILY ADJUST  
HANDLE HEIGHT**

**Code**  
**091558360**

**Description**  
Breeze Indoor Rollator

**Colour**  
○ White

100 kg  
15¾ st 220½ lb  
Maximum  
user weight

**1** YEAR  
Product  
warranty

**2** YEAR  
Main frame  
warranty


# Upright Walker

## Revolutionary mobility with forearm support


The Days Upright Walker mobility solution was designed to promote safe and secure posture for users. Crafted from sleek black aluminium, the walker boasts a robust design capable of supporting a maximum user weight of 330 lb. Minimal assembly is needed to ensure user-friendly setup, while practical features include a convenient cane holder and a showerproof bag for transporting personal belongings. Safety is paramount with loop lockable brakes ensuring a secure walk.

## Key features

- Large detachable storage back for carrying belongings
- Wheels are made from TPR rubber
- Cane/walking stick holder



|                                                   |                                                                                                |
|---------------------------------------------------|------------------------------------------------------------------------------------------------|
| <b>A</b>                                          | Handle height<br>98-121 cm (38 <sup>3</sup> / <sub>4</sub> -47 <sup>5</sup> / <sub>8</sub> " ) |
| <b>B</b>                                          | Overall depth<br>88 cm (33.5")                                                                 |
| <b>C</b>                                          | Overall width<br>60 cm (23.5")                                                                 |
| <b>D</b>                                          | Wheel diameter<br>Front: 30.4 cm (12")<br>Rear: 20.3 cm (8")                                   |
| Seat size<br>45.7 (18") wide x<br>28cm (11") deep |                                                                                                |



**MOULDED ARM  
RESTS FOR EXTRA  
COMFORT**



**CANE/WALKING  
STICK  
HOLDER**



**GRIP HANDLES TO  
AID GETTING IN  
AND OUT OF THE  
SEAT**

**Code**  
**7028116**

**Description**  
Upright Walker

**Colour**  
● Black

136 kg  
21 1/2 st 300 lb  
Maximum  
user weight


1 YEAR  
Product  
warranty

# Steel Tri-wheel Walker

## Sturdy, robust and economical

Steel Tri-wheel Walkers are designed to increase stability while walking with a three-wheel style that improves manoeuvrability without compromising support. The loop lockable breaks with adjustable tension is an important safety feature, along with the ergonomic handgrips and height-adjustable handles for comfort.

- **Folds easily for transport and storage**
- **Available with push-down or lockable loop brakes**
- **Steel frame with chrome or painted finish**



|                                  |                                      |
|----------------------------------|--------------------------------------|
| <b>A</b>                         | Handle height<br>83-93 cm (32¾-36½") |
| <b>B</b>                         | Overall depth<br>61 cm (24")         |
| <b>C</b>                         | Overall width<br>68.5 cm (27")       |
| <b>D</b>                         | Wheel diameter<br>19 cm (7½")        |
| Product weight<br>7.7 kg (17 lb) |                                      |



| Code      | Description      |
|-----------|------------------|
| 091356757 | Tri Wheel Walker |
| 091325505 | Tri Wheel Walker |
| 091356765 | Tri Wheel Walker |
| 091439355 | Tri Wheel Walker |
| 091356724 | Basket and Tray  |

| Colour   |
|----------|
| ● Blue   |
| ● Chrome |
| ● Red    |
| ● Grey   |

125 kg

19½ st 27½ lb

Maximum user weight

1 YEAR

Product warranty

2 YEAR


Main frame warranty

# Aluminium Lightweight Tri-wheel Walker

## Light and easy to manoeuvre

Aluminium Lightweight Tri-wheel Walkers feature a lightweight aluminium frame that is easy to manoeuvre, making them ideal for use indoors or outdoors.

- Lockable loop brakes for security and convenience
- Available with or without a basket
- Includes durable bag to transport and conceal your belongings



|                                  |                                      |
|----------------------------------|--------------------------------------|
| <b>A</b>                         | Handle height<br>83-93 cm (32¾-36½") |
| <b>B</b>                         | Overall depth<br>61 cm (24")         |
| <b>C</b>                         | Overall width<br>68.5 cm (27")       |
| <b>D</b>                         | Wheel diameter<br>19 cm (7½")        |
| Product weight<br>7.7 kg (17 lb) |                                      |



### With Bag

091553015  
091356732  
091553007  
091552728  
091541929  
091552710

Accessories  
091555242

### With Bag and Basket

091557461  
091557479  
091557487  
091557495  
091557503  
-

### Description

Aluminium Tri-Wheel Walker  
Aluminium Tri-Wheel Walker  
Aluminium Tri-Wheel Walker  
Aluminium Tri-Wheel Walker  
Aluminium Tri-Wheel Walker  
Aluminium Tri-Wheel Walker

Basket and Tray

### Colour

• Ruby Red  
• Blue  
• Quartz  
• Graphite  
• Purple  
• Russet Orange



OPTIONAL BASKET AND TRAY

100 kg  
15¾ st 220½ lb  
Maximum user weight

1 YEAR  
Product warranty

2 YEAR  
Main frame warranty



# Crutches

All Days Crutches are made from lightweight aluminium with a tough anodised finish, making them incredibly durable and easy to use.

They are designed for minimum maintenance and only require an occasional wipe-down with a soft cloth. Disinfectant solution can be used if needed.



# Crutches

---

## How to select crutches

Factors include handle options, adjustment levels and standard features

### Handle options

---



#### Standard

Standard plastic handgrips are easily wiped clean to reduce the risk of cross-contamination.



#### Ergonomic

Ergonomic handles were designed and developed with physiotherapists and healthcare professionals. The unique handle distributes pressure over a greater area of the user's hand, improving comfort and safety.



#### Comfy grip

Comfy grip handles are made from soft-touch plastic to provide a more comfortable and cushioned grip — ideal for long-term use.


### Adjustment levels

---



#### Single adjustable crutches

Single adjustable crutches enable you to adjust the handle height, but the distance between the handle and arm cuff is fixed.



#### Double adjustable crutches

Double adjustable crutches enable you to adjust the handle height and the distance between the handle and arm cuff.



#### Permanent crutches

Permanent crutches are not adjustable and must be cut to size to suit each user. They are ideal for long-term use.

### Standard features

---



#### Closed fixed armband

The closed fixed armband enables you to retain the crutch on your arm even if you release your grip on the handle. This prevents the crutch from dropping to the floor accidentally.



#### Non-marking ferrules

Non-marking ferrules are slip-resistant for added safety. The 22-mm ferrules are fitted for adult crutches, while 16-mm ferrules are fitted for children's crutches.



#### Push-pin adjustable

Adjustment is easy using the simple push-pin operation.

# Standard Grip Elbow Crutches\*

|          |                   |
|----------|-------------------|
| <b>A</b> | Handle to armband |
| <b>B</b> | Handle to ground  |



|                          |            |                             |                            |
|--------------------------|------------|-----------------------------|----------------------------|
| <b>Single adjustable</b> |            | <b>A: Handle to armband</b> | <b>B: Handle to ground</b> |
| <b>091275965</b>         | Standard   | 23 cm (9")                  | 66-94 cm (26-37")          |
| <b>091356468</b>         | Extra long | 23 cm (9")                  | 83-110 cm (32½-43½")       |
| <b>Double adjustable</b> |            | <b>A: Handle to armband</b> | <b>B: Handle to ground</b> |
| <b>091325323</b>         | Children's | 17-22 cm (6¾-8¾")           | 66-94 cm (26-37")          |
| <b>091115799</b>         | Standard   | 20-27.5 cm (8-10½")         | 66-94 cm (26-37")          |
| <b>091356443</b>         | Extra long | 23 cm (9")                  | 83-110 cm (32½-43½")       |

\*Crutches are sold as pairs.




# Ergonomic Grip Elbow Crutches\*

|          |                   |
|----------|-------------------|
| <b>A</b> | Handle to armband |
| <b>B</b> | Handle to ground  |



|                                   |            |                             |                            |
|-----------------------------------|------------|-----------------------------|----------------------------|
| <b>Double adjustable</b>          |            | <b>A: Handle to armband</b> | <b>B: Handle to ground</b> |
| <b>091356484</b>                  | Standard   | 20-27.5 cm (8-10½")         | 69-94 cm (27-37")          |
| <b>091356492</b>                  | Extra long | 20-27.5 cm (8-10½")         | 83-110 cm (32½-43½")       |
| <b>Permanent use, cut to size</b> |            | <b>A: Handle to armband</b> | <b>B: Handle to ground</b> |
| <b>091356476</b>                  | Standard   | 23.5-31 cm (9¼-12¼")        | 103 cm (40½")              |

\*Crutches are sold as pairs.



# Comfy Grip Elbow Crutches\*

|          |                   |
|----------|-------------------|
| <b>A</b> | Handle to armband |
| <b>B</b> | Handle to ground  |



Single adjustable  
**091439199** Standard  
**091531540** Extra long


Double adjustable  
**091439207** Standard  
**091531557** Extra long

**A: Handle to armband**  
 23 cm (9")  
 23 cm (9")

**A: Handle to armband**  
 66-94 cm (26-37")  
 83-110 cm (32½-43½")

**B: Handle to ground**  
 69-94 cm (27-37")  
 83-110 cm (32½-43½")

**B: Handle to ground**  
 20-27.5cm (8-10½")  
 20-27.5cm (8-10½")



\*Crutches are sold as pairs.

# Arthritis Crutches\*

|          |                  |
|----------|------------------|
| <b>A</b> | Handle to ground |
|----------|------------------|



Single adjustable  
**091437888**

**A: Handle to ground**  
 94-119 cm (37-47")



\*Crutches are sold as pairs.



# Walking Sticks

Select the walking stick that is right for you

Days Walking Sticks are available in various styles. Use the guide on the following page to help you select the right walking stick for you.




# Walking Sticks

## How to select walking sticks

Factors include stick height, left- or right-handed and type

### Walking stick height




Using a walking stick at the correct height is essential. When measuring, the handle should rest on the floor while your arm is relaxed. Then, adjust or trim the stick to the crease of the wrist.

When using your walking stick, you should maintain an upright posture with the elbow slightly bent. Your body weight will be supported when you push down on the stick while walking.

If you need two walking sticks, the length should be slightly longer as they will be used ahead of the body when you are walking.


### Left or right hand




You should hold the stick in hand on the opposite side to the one that needs support if you are using a walking stick to compensate for an injury to your leg or joint on one side.

If you use the stick for balance, you can use it with either hand.

### Types of walking sticks



**Adjustable**  
Adjustable walking sticks can be modified to different heights using simple spring-loaded catches and telescopic shafts.



**Tripod and tetrapod**  
Tripod and tetrapod walking sticks have a broad footprint to provide additional stability. These are best suited for heavier users who require a more stable base.



**Folding walking sticks**  
Folding walking sticks are lightweight metal sticks with sectioned shafts that enable you to fold them for storage. For example, you may transport them within a handbag. Strong elastic runs inside the shaft to ensure that the stick remains stable in its open position.


# Adjustable Aluminium Walking Stick

**A** Handle to ground

125 kg

19½ st 275½ lb

Maximum user weight



|                  |             |                            |
|------------------|-------------|----------------------------|
| <b>Code</b>      | <b>Type</b> | <b>A: Handle to ground</b> |
| <b>091565951</b> | Standard    | 71-96.5 cm (28-38")        |
| <b>091565969</b> | Long        | 84-109 cm (33-43")         |


# Standard Adjustable Walking Stick

**A** Handle to ground

125 kg

19½ st 275½ lb

Maximum user weight



|                  |             |                            |
|------------------|-------------|----------------------------|
| <b>Code</b>      | <b>Type</b> | <b>A: Handle to ground</b> |
| <b>091325588</b> | Standard    | 74-99 cm (29-39")          |

# Adjustable Aluminium Walking Stick

**A** Handle to ground

127 kg

20 st 280 lb

Maximum user weight



|                |             |                            |
|----------------|-------------|----------------------------|
| <b>Code</b>    | <b>Hand</b> | <b>A: Handle to ground</b> |
| <b>7500252</b> | Left        | 71-96.5 cm (28-38")        |
| <b>7500253</b> | Right       | 74-96.5 cm (29¼-38")       |


# Ergonomic Handle Walking Stick

**A** Handle to ground

125 kg

19½ st 275½ lb

Maximum user weight



|                  |             |                            |
|------------------|-------------|----------------------------|
| <b>Code</b>      | <b>Hand</b> | <b>A: Handle to ground</b> |
| <b>091325513</b> | Left        | 73-96 cm (28½-37½")        |
| <b>091325521</b> | Right       | 73-96 cm (28½-37½")        |

## Folding Steel Walking Stick

**A** Handle to ground

125 kg  
19½ st 275½ lb  
Maximum user weight




**Code** **A: Handle to ground**  
**091356914** 85-95 cm (33-37")  
**091356922** 74-83.5 cm (29-33")

## Quadruped Walking Stick

**A** Handle to ground

125 kg  
19½ st 275½ lb  
Maximum user weight



**Code** **Base Size A: Handle to ground**  
**091115344** Small 69-94 cm (27-37")  
**091356971** Wide 74-100 cm (29-39")

## Tripod Walking Stick

**A** Handle to ground

125 kg  
19½ st 275½ lb  
Maximum user weight



**Code** **Base Size A: Handle to ground**  
**091356930** Small 69-94 cm (27-37")  
**091325539** Wide 74-100 cm (29-39")

**A** Handle to ground

**B** Seat to ground

100 kg  
15¾ st 220½ lb  
Maximum user weight



**Code** **A: B:**  
**091115658** 84.5-95 cm (33¼-37½") 45.5-54 cm (18-21¼")



## Performance Health

### INTERNATIONAL SALES OFFICE

208 Kifissias Avenue, 15231 Halandri - Athens - Greece

tel.: +30 210 6849 970, +30 210 6712 266

fax: +30 210 6829 255

[info@performancehealth.com](mailto:info@performancehealth.com)

### CORPORATE HEADQUARTERS

1245 Home Avenue, Akron - Ohio 44310 - USA

tel.: +1 330 633 8460



**PERFORMANCE  
HEALTH™**

[performancehealth.global](http://performancehealth.global)

