



PIVOTAL HEALTH SOLUTIONS INTERNATIONAL PRODUCT GUIDE



TABLE OF CONTENTS

ABOUT US.....	Page	5
ERGOSTYLE™ ES2000 GEN II.....	Page	6
ERGOSTYLE™ FX.....	Page	8
ERGOSTYLE™ ELEVATING HYLO.....	Page	10
DOC DECOMPRESSION TABLE.....	Page	12
ERGOBASIC SERIES.....	Page	14
ERGOBASIC WITH MANUAL PUMP ELEVATION.....	Page	15
ERGOBENCH SERIES.....	Page	16
ATT-300 INTERSEGMENTAL ROLLER MASSAGE TABLE.....	Page	17
ERGOWAVE ROLLER MASSAGE TABLE.....	Page	18
THERA-P™ TRACTION TABLE.....	Page	20
ELEVATING THERAPY TILT TABLE.....	Page	22
LEG & SHOULDER (LAST) TABLE	Page	23

THERA-P™ ELECTRIC TREATMENT TABLE.....	Page 24
ESSENTIAL THERA-P™ ELECTRIC TREATMENT TABLE.....	Page 26
THERA-P™ BARIATRIC ELECTRIC TREATMENT TABLE	Page 27
EPD EPIDURAL POSITIONING DEVICE	Page 28
E-EPD ELECTRIC EPIDURAL POSITIONING DEVICE	Page 29
PT2000 BARIATRIC HI-LO MAT TABLE.....	Page 30
ALUMA ELITE BASIC TREATMENT TABLE.....	Page 31
PORTABLE ALUMINUM TREATMENT TABLE.....	Page 32
CLASSIC SERIES ATHENA MASSAGE TABLE.....	Page 33
STAIRCASES	Page 36
FULLY ELECTRIC PLATFORM MOUNT PARALLER BARS.....	Page 38
PORTABLE PARALLER BARS.....	Page 40
STOR-EDGE™ MULTI PURPOSE CARTS.....	Page 42
AMBULATION TRAINING.....	Page 46
CLINICAL ACCESSORIES.....	Page 47



ABOUT US

Pivotal Health Solutions is recognized for innovative design, exceptional craftsmanship, and ergonomically advanced healthcare equipment. Headquartered in Watertown, South Dakota, **United States**, our mission is to develop solutions that empower healthcare professionals to deliver better care while enhancing comfort, efficiency, and clinical performance.

From **chiropractors and physical therapists to athletic trainers, anesthesiologists, massage therapists, and occupational therapists**, healthcare professionals around the world rely on our extensive portfolio of products to support and grow their practices.

As the healthcare landscape advances, Pivotal Health Solutions remains committed to delivering reliable, high-performance products that support professionals and improve patient care worldwide.

Since its launch in 2009, **PHS Chiropractic** - a brand of Pivotal Health Solutions - has continued to evolve, developing some of the most comprehensive chiropractic adjusting table lines available today. Known for beautifully designed and highly customizable **Chiropractic tables** and accessories, the brand also offers the advanced **DOC Decompression System** designed to support a wide range of healthcare professionals.

PHS Medical designs and manufactures innovative therapy, medical, and pain management products built on the highest standards of safety, reliability, and quality. Our product range includes advanced physical therapy equipment, the patented **Epidural Positioning Device (EPD)**, the **Strava Splint Bath**, and a variety of high-quality **Parallel bars**. Each of these are designed to support healthcare professionals, improve patient outcomes, and add lasting value to medical facilities.


ERGOSTYLE™ ES2000 GEN II

The ErgoStyle™ ES2000 Gen II takes the table you know one step further. With smooth electric elevation and easy upgrades, you can customize your table to your practice setting. Our patented ErgoStyle drop elevates across the entire cushion, rather than a single pin, allowing for improved accuracy and a smoother adjustment with less force-for less stress on you and your patient. It's also fully modular, allowing for easy modifications and upgrades in the field, both now or later as your techniques and practice grow.

Standard Features:

- Lifting capacity: 181 kg
- Overall dimensions: 178cm x 71cm x 51-81cm
- Extended maximum length: 193cm with ankle extension
- Electric elevation from 51-81cm
- Single foot pedal for elevation
- Expanding covers over lifting mechanism
- Tilting headpiece with adjustable sewn face cushions
- 57cm Wide seamless cushions with multi-density foam
- Locking 15cm extending ankle rest
- Stationary chest
- Fixed pelvic section
- Paper roll holder and cutter
- Durable powder coat finish
- Dual auto-cocking foot pedal included with auto-cocking drops
- ADA compliant
- MDR compliant
- Rated voltage 120V/240V
- Rated frequency 50Hz/60Hz





ErgoStyle™ ES2000 GEN II
ES9300 shown with options

LIFTING CAPACITY
181 kg

OVERALL DIMENSIONS
178cm x 71cm x 51-81cm

EXTENDED MAXIMUM LENGTH
183cm with ankle extension

CUSHION WIDTH
57cm


ERGOSTYLE™ FX

The ErgoStyle™ FX offers full, ultra-smooth counterbalanced flexion and distraction, effortlessly meeting the core needs of your practice. Manual flexion/distraction is built right in, with infinitely variable electric spring tension so you can move the caudal section in any direction with just the touch of a finger. For added flexibility, the table also quietly elevates from its lowest height of 51cm to its full height of 81cm. It's also fully modular, allowing for easy modifications and upgrades in the field, both now or later as your techniques and practice grow.

Standard Features:

- Lifting capacity: 181 kg
- Overall dimensions: 178cm x 71cm x 51-81cm
- Extended maximum length: 198cm with ankle extension
- Electric elevation from 51-81cm
- Single foot pedal for elevation
- Expanding covers over lifting mechanism
- Tilting headpiece with adjustable sewn face cushions
- 57cm Wide seamless cushions with multi-density foam
- Manual flexion/distraction (Flexion: 15°; Extension: 12°; Side Bending: 30° total)
- Single tension foot pedal (manual flexion only)
- Locking 20cm extending ankle rest with ankle cuffs
- Stationary chest
- Paper roll holder and cutter
- Durable powder coat finish
- Single auto-cocking foot pedal included with auto-cocking drops
- ADA compliant
- MDR compliant
- Rated voltage 120V/240V
- Rated frequency 50Hz/60Hz





ErgoStyle™ FX
ES5820 shown with options

LIFTING CAPACITY
181kg

OVERALL DIMENSIONS
178cm x 71cm x 51-81cm

EXTENDED MAXIMUM LENGTH
198cm with ankle extension

CUSHION WIDTH
57cm

ERGOSTYLE™ ELEVATING HYLO

ErgoStyle™ Elevating HYLO is feature-rich with smooth elevation and up to 63° tilting capabilities that reduce stress on you and your patients. The ErgoStyle Elevating HYLO Table also features a tuck away foot plate allowing better access to the end of the table. Because it's fully modular, it accommodates easy modifications and upgrades as your practice and techniques grow. Choose from manual drops, various adjustable headpieces, breakaway chest, and many other options to customize your table to your practice needs.

Standard Features:

- Lifting capacity: 181kg
- Overall dimensions: 178cm x 71cm x 51-76cm
- Extended maximum length: 193cm with ankle extension
- Electric elevation from 51-76cm
- 0-63° Tilting
- Single foot pedal for elevation & HYLO
- Tilting headpiece with adjustable sewn face cushions
- 57cm Wide seamless cushions with multi-density foam
- Locking 15cm extending ankle rest
- Stationary chest
- Fixed pelvic section
- Foldaway foot plate
- Paper roll holder and cutter
- Durable powder coat finish
- Dual auto-cocking foot pedal included with auto-cocking drops
- ADA compliant
- MDR compliant
- Rated voltage 120V/240V
- Rated frequency 50Hz/60Hz





ErgoStyle™ Elevating HYLO EH9400
shown with options



LIFTING CAPACITY
181kg

OVERALL DIMENSIONS
178cm x 71cm x 51-76cm

EXTENDED MAXIMUM LENGTH
193cm with ankle extension

CUSHION WIDTH
57cm


DOC DECOMPRESSION TABLE

The DOC Decompression Table takes your practice to the next level with its unique design and impressive list of standard features. Suitable for physiotherapists, chiropractors and doctors. The state-of-the-art digital command center controls specific vertebral targeting, including axial rotation and lateral flexion and correctly positions patients in antalgic posture. Offers separate lumbar and cervical decompression programming and delivers continuous readout and graphing of treatment protocols. The digital command center also controls table elevation and rear flexion, all at the push of a button. Patients may be treated in supine or prone positions and from various angles.

Standard Features:

- Lifting capacity: 159kg
- Overall dimensions: 76cm x 196cm x 53cm
- Electric elevation range: 53-74cm
- Lateral flexion positioning and axial rotation
- Cervical capture headpiece
- Traction belting system enhances restraint for optimal decompression
- User-defined treatment hold and relax times
- Real-time digital treatment tracking
- Speed control
- Pre-programmed lumbar and cervical decompression protocols
- Digital Command Center includes: Specific vertebral targeting, separate lumbar and cervical traction/decompression programming, continuous readout and graphing of treatment parameters, and interactive color display
- Reference library - including patient setup and training materials
- MDR compliant
- Rated voltage 120V/240V
- Rated frequency 50Hz/60Hz





Supine Flexion and Extension

The ideal position for older or heavier patients who can't tolerate decompression but benefit from extension with disc herniations or disc bulging.



Prone Flexion and Extension

The ideal patient position for younger patients who can tolerate prone positioning for symptoms of disc herniations and disc bulges.



Lateral Flexion

Laterally flexing away from a disc herniation can dramatically help improve disc herniation outcomes. It is also effective for symptoms of antalgic patients.



Rotation

Effective patient positioning for antalgic patients.



E9011-B

- FDA Cleared
- EU Certified
- Licensed in Canada

LIFTING CAPACITY
159kg

OVERALL DIMENSIONS
76cm x 196cm x 53cm

ELECTRIC ELEVATION RANGE
53-74cm

CUSHION WIDTH
76cm

ERGOBASIC SERIES

Standard Features:

- Load capacity: 159kg
- Overall dimensions: 69cm x 165cm x 46-61cm
- Tilt headpiece
- Accelerator™ III pelvic drop
- Unitized steel frame
- 52cm Wide seamless cushions
- Paper roll holder and cutter
- Shipped fully assembled
- ADA Compliant
- MDR Compliant

The ErgoBasic Series offers a critical mix of adjusting features that is praised by many practitioners.

ErgoBasic.....EB9040
(Soft Foam)

ErgoBasic.....EB9050
(Firm Foam)



EB9050



EB9050 with options

ERGOBASIC WITH MANUAL PUMP ELEVATION

Standard Features:

- Load capacity: 181kg
- Overall dimensions: 69cm x 177cm x 46-61cm or 56-71cm
- Manual hydraulic pump elevation
- Tilting headpiece with adjustable sewn face cushions
- Firm foam
- Paper roll holder and cutter
- Durable powder coat finish
- Accelerator™ III Pelvic Drop
- ADA Compliant
- MDR Compliant

ErgoBasic.....EB9060
(46-61cm)

ErgoBasic.....EB9061
(56-71cm)



EB9060

ERGOBENCH SERIES

Standard Features:

- Load capacity: 159kg
- Overall dimensions: 61cm x 185cm x 46-61cm
- Unitized steel frame
- 61cm Seamless cushions
- Multi-fit arm rests
- Paper roll holder and cutter
- Bottleneck shape
- ADA Compliant
- MDR Compliant

The ErgoBench Series brings uncompromised stability and style to any office at an affordable price. It features a smooth, flat surface, unitized steel frame with a bottleneck shape, seamless cushions, multi-fit arm rests, and a paper roll holder and cutter.

Fixed Top...EB9020-U
(Soft Foam)

Tilt Head...EB9000-U
(Soft Foam)



EB9020



EB9000

ATT-300 INTERSEGMENTAL ROLLER MASSAGE TABLE

Standard Features:

- Load capacity: 156kg
- Overall dimensions: 56cm x 183cm x 64cm
- Clockwise and counter-clockwise rotation
- Travel indicator
- Automatic lowering device
- Adjustable massage pressure
- Vibration
- LED roller height indicator
- 30 Minute timer
- (1) Cervical bolster and (1) knee bolster (black)
- 10cm Roller height adjustment
- Four ball casters
- Rated voltage 120V

ATT-300-W

This Wood Roller Massage Table features three large rotating rollers traveling both clockwise and counterclockwise, up and down the spine, to knead and relax muscles in the neck, shoulders and lower back.



ATT-300



ATT-300-W

ERGOWAVE ROLLER MASSAGE TABLE

Standard Features:

- Load capacity: 159kg
- Overall dimensions: 198cm x 64cm x 61cm
- Quick height adjustments: Gas spring makes roller height adjustments quick and easy
- Ultra quiet operation: high-quality motor for nearly silent operation without disturbing vibrations
- Built-in storage shelf
- Welded steel frame for stability
- 30-Minute timer
- 8cm Roller height adjustment

The ErgoWave Roller Massage Table covers all areas of the spine isolating the lumbar, thoracic or cervical regions. Three independent rollers repeat every 13 seconds delivering a unique wave action that travels up the spine. A powerful yet nearly silent motor never disturbs your patient's relaxation and you can add a leg bolster or head pillow for added comfort.

ErgoWave Roller Massage.....EW9080



EW9080

CUSTOMIZE IT!

YOU CAN CHOOSE COLORS AND BRANDING FOR MOST PRODUCTS THAT ARE DESIGNED TO MATCH FOR A MORE COHESIVE LOOK IN YOUR FACILITY.



THERA-P™ TRACTION TABLE

Comprised of 4 or 6 sections, the Thera-P Traction Tables use electric power to help with manual adjustment, allowing for ease of use and control of table positioning. For the practitioner, the Thera-P Traction Table offers a 51-94cm full height adjustment and section mobility so the patient can be placed in the proper position for maximum results with minimum effort. Each section is adjustable, providing numerous achievable patient positioning options.

Standard Features:

- Load capacity: 294kg
- 69cm x 198cm x 51-94cm
- Durable and stylish steel base
- White powder coat base
- Whisper quiet lift system
- Height adjustment: 51-94cm
- Dual foot bar control
- 69cm Wide cushions in 2 different top configuration options
- 5cm High density foam cushions
- Rounded corners
- Friction-free gas shock cushion separation with lock and unlock lever
- (2) Low profile locking casters
- Head section angle range: +30° to -90°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 18cm
- Lift back: 0 to 70°
- MDR compliant



Knee & Cervical Bolster Set
and Traction Belt



Traction Head



Cervical Attachment
& Head Pillow



TXT-1004 with options

**Branding Options
Available**



TXT-1004



TXT-1006

ELEVATING THERAPY TILT TABLE

The Elevating Therapy Tilt Table is ideal for nursing home and rehabilitation facility use. Smooth and secure elevation adjusts from 57-100cm high, and tilts from a horizontal position to 90°. A heavy-duty powder coat frame supports a load capacity of 159kg and includes an extra-large foldaway foot plate with non-skid surface and dual safety straps for added stability. The 5cm high density foam cushions with two configurable options are designed for patient comfort.

Standard Features:

- Load capacity: 159kg
- 69cm x 208cm x 57-100cm
- (2) Safety straps
- Extra large foldaway foot plate with non-skid tape
- Hand control for tilt
- Single foot control for lift
- Nose hole with plug
- (4) Individual locking casters



LEG & SHOULDER THERAPY (LAST) TABLE

Standard Features:

- Load capacity: 294kg
- 76cm x 213cm x 48-94cm
- Adjustable top with standard elevating base
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system
- Height adjustment: 48-94cm
- Single foot control pedal
- Elevating back: 0 - 70°
- Elevating/drop legs: 0 - 90°
- Leg height: 43cm
- 5cm High density foam cushions
- Nose hole
- Rounded corners
- Low profile locking casters
- ADA Compliant
- MDR compliant



PT9060 with options
Patent # 10,857,052

THERA-P™ ELECTRIC TREATMENT TABLE

Ease of use, multiple options and a sleek but durable design make the THERA-P line of electric treatment tables the choice of therapists everywhere. Seven top options allow therapists to choose the table that best meets their needs and that ensures patients' most comfortable positioning. The THERA-P electric treatment table is constructed of precision-cut flat and tube steel for durability and style. It also has a whisper-quiet lift system (with elevating foot switch) and easy-to-operate caster system for better ergonomics.

Standard Features:

- Load capacity: 295kg
- 69cm x 196cm x 43-89cm
- Durable and stylish steel base
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 43-89cm
- 69cm Wide cushions in 7 different configuration options
- 5cm High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90° (PT1177: +35° to -45°)
- Lower extremity angle range: 0 to 65° (PT1177: 0-75°)
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 18cm
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant
- MDR compliant
- Rated voltage 120V/240V
- Rated frequency 50Hz/60Hz



1 Section Thera-P
PT1101/EPT01



2 Section Thera-P
PT1102/EPT02



3 Section Thera-P
PT1103/EPT03



4 Section Thera-P
PT1104/EPT04



5 Section Thera-P
PT1105/EPT05



6 Section Thera-P
PT1106/EPT06



7 Section Thera-P
PT1177/EPT77

Options:

240V Tables available

PT303x - 76cm Cushion Configuration

S1135 - Manual Elevating Midsection (Only available on EPT03 & EPT05)

S1077 - Dual Foot Control Pedal

S1078 - Dual Foot Bar Control

S2078 - Dual Foot Bar Control with Front Foot Bar Adder

EPTPH01 - Paper Roll Holder

(Only available on EPT01, EPT02, or EPT03)

PT10001 - Nose Hole Plug

S1151 - (1) Battery and (1) wall charger

S1152 - (2) Batteries and (1) wall charger



ESSENTIAL THERA-P ELECTRIC TREATMENT TABLE

Standard Features:

- Load capacity: 295kg
- 69cm x 196cm x 43-89cm
- Durable and stylish steel base
- White anti-bacterial and anti-microbial powder coat base with single foot control pedal
- Whisper quiet lift system
- Height adjustment: 43-89cm
- 69cm Wide cushions in 7 different configuration options
- 5cm High density foam cushions with rounded corners
- Individual locking casters
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant
- MDR compliant
- Rated voltage 120V/240V
- Rated frequency 50Hz/60Hz



EPT03

THERA-P™ BARIATRIC ELECTRIC TREATMENT TABLE

Standard Features:

- Load capacity: 295kg
- 86cm or 91cm x 196cm x 43-89cm
- Durable and stylish steel base
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 86cm or 91cm Wide cushions in 4 different configuration options
- 5cm High density foam cushions with rounded corners
- Low profile locking casters
- ADA Compliant



PT3613 with options

EPD EPIDURAL POSITIONING DEVICE

Patented in 2001, our EPD is a landmark device designed by an anesthesiologist that positions patients correctly and comfortably by encouraging cervical, thoracic and lumbar flexion while maintaining a solid and stable position. Manually controlled chest, arm and head support ensures that patients of every body type and size can be accommodated and the back remains immobile during epidural placement. The EPD is suitable for use in operating rooms, labor and delivery suites, pain centers or anywhere epidural and spinals are performed.


Standard Features:

- Load capacity: 272kg
- 52cm x 76cm x 127-152cm
- White powder coat base
- 180° Adjustable face rest
- 6 Position adjustable arm rest
- Arm rest size: 42cm x 25cm
- Arm rest adjustment: 81-109cm from the floor
- Interchangeable torso cushions (1 flat and 1 elevated)
- Torso cushion size: 20cm x 20cm
- Torso cushion adjustment: 14-25cm
- MDR compliant



FRC 30-50

EPD
(Epidural Positioning Device)



E-EPD ELECTRIC EPIDURAL POSITIONING DEVICE

Based on a design by an anesthesiologist, the E-EPD positions patients correctly and comfortably. A new curved design accommodates patients of a variety of body sizes and weights. Electronic controls adjust the E-EPD to the correct height and foot position for each patient, while the easy mechanical lock and release system makes chest, arm and head adjustments quick and effortless.

Standard Features:

- Load capacity: 272kg
- 71cm x 74cm x 127-152cm
- 180° Adjustable face rest
- Prepositional pillow for use during cervical, lumbar and thoracic placement
- Arm rest adjustment: 15cm
- Torso cushion adjustment: 15cm
- Single hand control adjusts both midsection and foot rest
- (2) Rechargeable batteries last up to 100 cycles each
- Medical grade profile locking casters
- Rated voltage 120V



E-EPD
(Electric Epidural
Positioning Device)

PT2000 BARIATRIC HI-LO MAT TABLE

The PT2000 Bariatric Hi-Lo Mat Table is ideal for the rehabilitation setting and is available in a range of sizes to provide your facility with the option of delivering optimum care for larger patients.

Standard Features:

- Load capacity: 454kg
- Available in 1.22m x 2.13m, 1.52m x 2.13m, 1.83m x 2.44m configurations
- Height range: 46-84cm
- Durable steel base with powder coat finish
- 4cm High density foam cushion
- Rounded corners
- Single foot control
- Convenient drop down casters
- ADA Compliant



ALUMA ELITE BASIC TREATMENT TABLE

The dependable Aluma Elite Basic Treatment Table is the perfect centerpiece for each therapy or rehab facility. The extraordinary load capacity of 295kg and the welded aluminum frame make the Aluma Elite Basic Treatment Table the perfect table for your office, clinic, gym, or facility. Make it your own by adding length or width, laminate shelves or a cabinet that can be mounted on either end.

Standard Features:

- Load capacity: 295kg
- 69cm or 76cm x 183cm or 198cm x 79cm
- Welded aluminum frame
- H-brace construction
- 5cm High density foam cushion
- Rounded corners
- Tough powder coat finish
- Available in Flat and Lift Back Configurations



A4327LB with options

PORTABLE ALUMINUM TREATMENT TABLE

The Portable Aluminum Treatment Table is rugged yet lightweight and functional. It comes standard with a 5cm firm foam system, rounded corners, adjustable height leg system, four field feet, face rest and standard carry case.

Standard Features:

- Load capacity: 272kg
- Table weight: 17kg
- Lightweight aluminum frame
- Height range: 61-84cm
- 11-Layer plywood base
- Auto-hold hinge
- 20cm on center dual end accessory ports
- 5cm Firm foam cushion, black or navy upholstery
- Rounded corners
- Removable field feet
- Black carry case
- Face rest and pillow



CLASSIC SERIES ATHENA MASSAGE TABLE

Standard Features:

- Load capacity: 204kg
- 71cm or 76cm x 185cm x 64-89cm
- Table weight: 14kg for 71cm, 15kg for 76cm
- Shiatsu release cables
- 23cm on center dual-end accessory ports Reiki end panels on both ends
- 8cm Multi-layer foam

Upgrades:

- True Touch Upholstery
- True Touch Upholstery with Touch Response Memory Foam
- Rounded Corners
- Breast Recesses
- Prenatal Upgrade: Breast Recesses and Belly Cutout with Sling




CLASSIC SERIES UPHOLSTERY CHOICES:



TRUE TOUCH UPHOLSTERY CHOICES:







ALUMINUM STAIRCASES



PTST-002-AHR-A-
Straight

Straight Staircases

Standard Features:

- Load capacity: 272kg
- Platform: 76cm x 61cm x 61cm
- (3) 76cm x 25cm x 15cm Steps
- Overall: 89cm x 155cm x 152cm
- 2.44m Ceiling height required
- Ships fully assembled



PTST-001-AHR-A-
Small Pass Through

Small Pass Through

Standard Features:

- Load capacity: 272kg
- Platform: 76cm x 61cm x 30cm
- (2) 76cm x 25cm x 10cm Steps
- (1) 76cm x 25cm x 15cm Steps
- Overall: 89cm x 155cm x 121cm
- 2.44m Ceiling height required
- Ships fully assembled



PTST-003-AHR-A -
Large Pass Through

Large Pass Through

Standard Features:

- Load capacity: 272kg
- Platform: 91cm x 91cm x 61cm
- (5) 91cm x 25cm x 10cm Steps
- (3) 91cm x 25cm x 15cm Steps
- Overall: 104cm x 314cm x 152cm
- 2.44m Ceiling height required

WOODEN STAIRCASES



PTST-005-WHR
corner configuration

Convertible Staircases

Standard Features:

- Load capacity: 159kg (Maple Laminate) & 204kg (Birch)
- Platform: 76cm x 76cm x 61cm
- (4) 76cm x 23cm x 15cm steps and (6) 76cm x 23cm x 10cm steps
- 2.74m Ceiling height required



PTST-001-WHR-ADA

Pass Through Staircases

Standard Features:

- Load capacity: 159kg (Maple Laminate) & 204kg (Birch)
- Platform is 76cm x 61cm x 30cm
- (2) 76cm x 23cm x 15cm steps and (3) 76cm x 23cm x 10cm steps
- 2.44m Ceiling height required



PTST-002-WHR

Straight Staircases

Standard Features:


- Load capacity: 159kg (Maple Laminate) & 204kg (Birch)
- Platform is 76cm x 61cm x 61cm
- (4) 76cm x 23cm x 15cm steps
- 2.74m Ceiling height required

FULLY ELECTRIC PLATFORM MOUNT PARALLEL BARS

Engineered for durability and precision, our Adjustable Parallel Bars support up to 362kg (188kg per bar) and are available in 2.4m, 3.1m, and 3.7m lengths. The heavy-duty steel platform features an anti-bacterial, anti-microbial powder coat finish for enhanced hygiene, while brushed stainless steel handrails ensure a comfortable grip.

The patented height and width adjustment system, powered by a microprocessor, offers independent or synchronized movement for ultimate flexibility. A user-friendly touch screen control panel displays bar position, ensuring safe and reliable operation. Designed for stability and safety, the platform is slip-resistant, making it ideal for various rehabilitation and training needs.





Standard Features:

- Load capacity: 362,8Kg (188Kg per bar)
- 2.4m, 3.1m, and 3.7m Lengths available
- Platform: 120cm W (48.6cm - 77.5cm W between columns)
- Heavy duty, rugged, patented height and width adjustment system
- Steel, multi-section platform with gray anti-bacterial and anti-microbial powder coat base
- Brushed stainless steel hand rails (5.1cm Diameter)
- Battery powered
- Digital touch screen display on conveniently located, user-friendly control panel shows exact bar height
- Memory buttons allow storage and recall of multiple bar height settings
- Microprocessor control unit ensures synchronous operation of all four columns
- Height monitoring functions ensure reliable operation and safe switch-off in case of a fall
- Bars can move independently or in unison
- Slip-resistant platform

PORTABLE PARALLEL BARS (PTPB1000)

Introducing a revolutionary, first-of-its-kind solution: Portable Parallel Bars designed for ultimate versatility and unmatched durability. Engineered with a 272.2kg base capacity and foldable design, they offer effortless transport and a compact footprint for storage. With adjustable aluminum handrails, dual-lock height adjustment, and a 55-degree turn radius, these bars adapt seamlessly to any environment. Built with a welded frame, anti-slip surface, and medical-grade casters, they're coated with anti-microbial protection to elevate safety and performance. Elevate rehabilitation and training with the innovation you've never seen before.

Standard Features:

- Load capacity: 272.2 Kg on the base, 136.1 Kg per bar
- Foldable with casters for easy transport
- Transport mode: 31.1 cm - 44.5 cm W x 152.4 cm L x 100.3 cm H
55° Turn radius
- Extended mode: 116.6 cm W x 152.4 cm L x 72.4 cm - 104.1cm H
- Adjustable handrail: 57.8 cm - 67.3 cm W
- Dual lock mechanism for height adjustment
- (4) 10.2 cm Total locking medical grade caster system
- Aluminum handrails
- Welded aluminum frame and walking surface
- Anti-microbial and anti-bacterial powder coat
- Anti-slip platform



STOR-EDGE™ MULTI PURPOSE CARTS



SR-001

Standard Features:

- 76cm x 46cm x 183cm
- Front: slatwall & band dispensers
- Back: slatwall
- Sides: weight racks
- (9) 10cm Double hooks
- (10) J-Hooks

Options:

- **SR-WK1** - Weight Kit 1
(2) .45-4.54kg Dumbbell sets
(1) 16 pc. Hugger weight set
(1) 5 Rep. resistance band set
- **SR-MI-001** - Mirror
50cm x 166cm (back side)



SR-002

Standard Features:

- 76cm x 61cm x 155cm
- Front: doors with (3) shelves
- Back: slatwall
- Sides: weight racks
- (4) 10cm Double hooks
- (10) J-Hooks

Options:

- **SR-WK2** - Weight Kit 2
(2) .45-4.54kg Dumbbell sets
(1) 16 pc. Hugger weight set
- **TSFO-FHLD**
Flush mount handles on large doors



Front

Back

SR-004

Standard Features:

- 61cm x 61cm x 155cm
- Front: band dispensers
- Back: slatwall
- Sides: weight racks
- (3) 10cm Double hooks
- (6) J-Hooks

Options:

- SR-WK4 - Weight Kit 4
 - (2) .45-4.54kg Dumbbell sets
 - (2) .45-2.27kg Dumbbell sets
 - (1) 5 Rep. resistance band set



Front

Back

SR-005

Standard Features:

- 32"W x 32"D x 72"H
- Front: cubbies, drop down desk with lock, (3) 20cm drawers
- Back: band dispensers & slatwall
- Left side: weight racks
- Right side: mirror (4) 10cm Double hooks (12) J-Hooks

Options:

- SR-WK3 - Weight Kit 3
 - (2) .45-4.54kg Dumbbell sets
 - (1) 16 pc. Hugger weight set
 - (1) 5 Rep. resistance band set
 - (1) .45-1.36cm Weight bar set
- TSFO-FHD
 - Flush mount handles on doors
- TSFO-FHDR
 - Flush mount handles on drawers

STOR-EDGE™ MULTI PURPOSE CARTS



Front Back

SR-006

Standard Features:

- 91cm x 91cm x 185cm
- Front: Doors with (2) shelves & (2) 25cm Drawers
- Flush mount handles
- Back: weight racks
- Sides: slatwall
- (4) 10cm Double hooks
- (12) J-Hooks

Options:

- **SR-WK3** - Weight Kit 3
 - (2) .45-4.54kg Dumbbell sets
 - (1) 16 pc. Hugger weight set
 - (1) 5 Rep. resistance band set
 - (1) .45-1.36kg Weight bar set



Front Back

SR-007

Standard Features:

- 76cm x 46cm x 152cm
- Front: slatwall
- Back: weight racks
- Right side: band dispensers
- Left side: storage compartment
- (4) 10cm Double hooks
- (10) J-Hooks

Options:

- **SR-WK3** - Weight Kit 3
 - (2) .45-4.54kg Dumbbell sets
 - (1) 16 pc. Hugger weight set
 - (1) 5 Rep. resistance band set
 - (1) .45-1.36kg Weight bar set



Front

Back

SR-008

Standard Features:

- 117cm x 66cm x 109cm
- Front: (5) 46cm drawers & (1) 46cm door
- Back: slatwall
- Left side: weight racks
- Right side: storage compartment
- (4) 10cm Double hooks
- (10) J-Hooks

Options:

- **TSFO-FHD**
Flush mount handles on doors
- **TSFO-FHLD**
Flush mount handles on large doors



Front

Back

SR-009

Standard Features:

- 91cm x 91cm x 185cm
- Front: doors with (2) Shelves & (2) 25cm drawers
- Back: weight racks
- Sides: slatwall
- (4) 10cm Double hooks
- (10) J-Hooks

Options:

- **SR-WK3** - Weight Kit 3
2) .45-4.54kg Dumbbell sets
(1) 16 pc. Huger weight set
(1) 5 Rep. resistance band set
(1) .45-1.36kg Weight bar set
- **SR-MI-009** - Mirror - 67.3cm x 166.4 (back side)
- **TSFO-FHD**
Flush mount handles on doors
- **TSFO-FHLD**
Flush mount handles on large doors

AMBULATION TRAINING



WCCT



NFS-100



WCRS

Ramp and Curb Training

Standard Features:

- Load capacity: 159kg
- 213cm x 112cm x 20cm
- Hardwood plywood construction
- Adjustable curbs store beneath platform
- Anti-slip safety tread on platform and all steps for increased safety

Nested Foot Stools

Standard Features:

- Load capacity: 159kg
- Overall: 41cm x 51cm x 20cm
- Hardwood plywood construction
- Anti-slip safety tread on platform for increased safety

Rolling Step

Standard Features:

- Load capacity: 136kg
- 76cm x 76cm x 20cm
- Hardwood plywood construction
- Anti-slip safety tread on platform for increased safety

CLINICAL ACCESSORIES



Smart Step Wellness Mats

Our line of high-performance Smart Step Wellness Mats deliver the comfortable support you need to help reduce stress, strain and fatigue while you work, play or exercise. Each mat features extra-durable construction, non-slip tops and bottoms, and an easy-to-clean, **antimicrobial** surface.

MDR Compliant



Rolling Stools

Our rolling stools combine comfort with control thanks to a 8cm deluxe seat cushion, a five-star base with heavy-duty casters and pneumatic controls that adjust seat height from 46-58cm or 51-66cm

MDR Compliant

Available options: RS, MRS, I9332, I9310



OMNI Cervical Relief Pillow

This pillow uses an ergonomically contoured shape to comfortably "cradle" the head and neck, allowing total relaxation of neck, shoulder and back muscles.

CEBL-1001
Omni Cervical Relief Pillow

PIVOTAL HEALTH SOLUTIONS INTERNATIONAL PRODUCT GUIDE

Performance Health

INTERNATIONAL SALES OFFICE

208 Kifissias Avenue, 15231 Halandri - Athens - Greece

tel.: +30 210 6849 970, +30 210 6712 266

fax: +30 210 6829 255

info@performancehealth.com

CORPORATE HEADQUARTERS

1245 Home Avenue, Akron - Ohio 44310 - USA

tel.: +1 330 633 8460

